



NORMAND INFO
Medical Solutions

REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

June 7, 2017

Over the past few days, many organizations around the world have fallen victim to the malicious “WannaCrypt” exploit.

Normand-Info provides a “Remisol Advance System Update” kit on a quarterly basis; however, the solutions for the vulnerability described in the Microsoft Security Bulletin MS17-010 are not included in the latest “Q1 2017” edition, and therefore require immediate action (instructions on the next page).

As a longer term solution, the “Remisol Advance System Update” kit will be upgraded in the following days, and will provide the complete set of patches to secure your system; we shall let you know in a separate communication letter.

Sincerely yours,

Jean-Louis GAVEL.

OFFICIAL MICROSOFT INFORMATION

The official information provided by Microsoft can be found at these links:

- Microsoft Security Bulletin MS17-010:

<https://technet.microsoft.com/en-us/library/security/ms17-010.aspx#Affected%20Software%20and%20Vulnerability%20Severity%20Ratings>

- Customer Guidance for WannaCrypt attacks:

https://blogs.technet.microsoft.com/msrc/2017/05/12/customer-guidance-for-wannacrypt-attacks/?utm_source=Direct



NORMAND-INFO

S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>

N° de SIRET 529 094 476 - Code APE 6201Z - RCS ARRAS 529 094 476 - N° Intra Communautaire FR 78529094476



REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

PROCEDURE

Please ensure you follow this procedure step-by-step!

Pre-requisite: these patches should be installed on up-to-date Remisol machines.

IDENTIFYING THE OPERATING SYSTEM (OS) OF THE REMISOL SERVER AND CLIENTS:

To find out what OS is currently running on the machine:

Start menu > Settings > Control panel > System and Security > System

Then please download the respective patch / Knowledge base for the OS – please read on.

DOWNLOADING THE MICROSOFT OPERATING SYSTEM (OS) PATCHES / KNOWLEDGE BASE (KB)

In regards to the Remisol Advance machines, whether hardware or virtualized, the following patches should be immediately installed:

OS	Patch / KB	File size
Windows XP	KB4012598	668 KB
Windows Server 2003	KB4012598	952 KB

- For Windows XP machines please download
✓ KB4012598:

http://download.windowsupdate.com/d/csa/csa/secu/2017/02/windowsxp-kb4012598-x86-custom-enu_eceb7d5023bbb23c0dc633e46b9c2f14fa6ee9dd.exe,

- ✓ Service Pack 3 package (prerequisite to KB install):

http://www.download.windowsupdate.com/msdownload/update/software/dflt/2008/04/windowsxp-kb936929-sp3-x86-enu_c81472f7eeee2eca421e116cd4c03e2300ebfde4.exe

- ✓ Executable (SYSADM.EXE prerequisite to KB install)

<http://www.remisoladvance.com/patches/sysadm.exe>



NORMAND-INFO

S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>



NORMAND INFO
Medical Solutions

REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

- For Windows Server 2003 machines please download:
 - ✓ KB4012598:

http://download.windowsupdate.com/c/csa/csa/secu/2017/02/windowsserver2003-kb4012598-x86-custom-enu_f617caf6e7ee6f43abe4b386cb1d26b3318693cf.exe

- ✓ Executable (SYSADM.EXE prerequisite to KB install)

<http://www.remisoladvance.com/patches/sysadm.exe>



NORMAND-INFO

S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>

N° de SIRET 529 094 476 - Code APE 6201Z - RCS ARRAS 529 094 476 - N° Intra Communautaire FR 78529094476



NORMAND INFO
Medical Solutions

REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

PREPARING A USB STICK WITH THE MICROSOFT OS PATCHES, SERVICE PACK AND SYSADM.EXE EXECUTABLE

To update Remisol machines, the OS patches will need to be downloaded to a USB flash drive.

On a machine connected to the Internet:

- Use a new USB stick with sufficient space for the OS patches needed. If downloading more than one OS patch, create a separate folder for each OS.
- Plug the USB stick to the machine.
- Scan the USB stick for viruses - **Please ensure your antivirus is up to date! If not please update it first before you scan the USB stick,**
- Connect to the Internet via your web browser,
- Download the adequate patch / KB /Service pack and the SYSADM.EXE executable for the Remisol machine OS that you are going to secure and save them onto the USB stick in their respective folder.
- Your USB stick is now ready.



NORMAND-INFO

S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>

N° de SIRET 529 094 476 - Code APE 6201Z - RCS ARRAS 529 094 476 - N° Intra Communautaire FR 78529094476



NORMAND INFO
Medical Solutions

REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

APPLYING THE MICROSOFT OS PATCHES

The Remisol production application will need to be stopped for the OS updates. Please coordinate with the lab manager with the understanding that this is a priority security update.

- **On a network of Remisol machines**

Here we mean Remisol machines connected to the hospital / laboratory network; for machines NOT connected to any network, please refer to **On standalone Remisol machines** sections below

Power off each of the Remisol client machines,

On the “Windows 2003 Server” Remisol server

On the Remisol server, stop all running applications and exit Remisol,

Log on the machine with the Remisol login with the highest administrative rights in your laboratory (sysadm)!

Back up the Remisol database. To do so, click on Windows **Start** button in the desktop lower left corner > (**All Programs** if applicable) > **3 – Save Remisol Advance database**; the Remisol database will be backed up to its default location: **D:\MSSQL\Backup** under the name **Custom.dat**.

Plug in the USB stick previously created,

Explore the USB stick folders and look for the one in which you have downloaded KB/Service Pack and the executable.

Double click on the file SYSADM.EXE, the Windows session will close.

Log on the machine again with the login with the highest administrative rights in your laboratory (sysadm)!

Double click on the patch / KB to launch its installation – follow onscreen instructions,

When completed, unplug the USB stick,

Restart the server.

The server is now patched.



NORMAND-INFO

S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>

N° de SIRET 529 094 476 - Code APE 6201Z - RCS ARRAS 529 094 476 - N° Intra Communautaire FR 78529094476



REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

On the Windows XP Remisol client(s)

Proceed with the Remisol client machines one by one.

- Restart the Remisol client,
 - Log on the machine with the login with the highest administrative rights in your laboratory (sysadm)!**
 - Do not start any application; in case Remisol automatically starts, stop it,
 - Plug in the USB stick,
 - Explore the USB stick folders and look for the one in which you have downloaded KB/Service Pack and the executable.
 - Double click on the file SYSADM.EXE, the Windows session will close.
 - Log on the machine again with the login with the highest administrative rights in your laboratory (sysadm)!**
 - Double click on the Service Pack 3 installer located in the USB stick and follow onscreen instructions.
 - Restart the Remisol client,
 - Log on the machine again with the login with the highest administrative rights in your laboratory (sysadm)!**
 - Double click on the patch / KB located in the USB stick to launch its installation – follow onscreen instructions,
 - When completed, unplug the USB stick,
 - Restart the client.
- The client is now patched.
- Repeat these operations for every client on the network.



NORMAND-INFO
S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>

Once all clients were patched, your network should now be immune to the WannaCrypt exploit.



NORMAND INFO
Medical Solutions

REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

- **On a Windows Server 2003 standalone Remisol machines i.e. not connected to the hospital / laboratory network**
 - Restart the standalone server,
 - Log on the machine with the Remisol login with the highest administrative rights in your laboratory (sysadm)!**
 - Back up the Remisol database. To do so, click on Windows **Start** button in the desktop lower left corner > (**All Programs** if applicable) > **3 – Save Remisol Advance database**; the Remisol database will be backed up to its default location: **D:\MSSQL\Backup** under the name **Custom.dat**,
 - Plug in the USB stick previously created,
 - Explore the USB stick folders and look for the one in which you have downloaded KB/Service Pack and the executable.
 - Double click on the file SYSADM.EXE, the Windows session will close.
 - Log on the machine again with the login with the highest administrative rights in your laboratory (sysadm)!**
 - Double click on the patch / KB to launch its installation – follow onscreen instructions,
 - When completed, unplug the USB stick,
 - Restart the standalone machine.

The standalone machine is now patched.



NORMAND-INFO

S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>

N° de SIRET 529 094 476 - Code APE 6201Z - RCS ARRAS 529 094 476 - N° Intra Communautaire FR 78529094476



NORMAND INFO
Medical Solutions

REMISOL ADVANCE SECURITY UPDATE FOR WIN XP/ WIN 2003

On a Windows XP standalone Remisol machines i.e. not connected to the hospital / laboratory network

- Restart the standalone Remisol client,
- Log on the machine with the login with the highest administrative rights in your laboratory (sysadm)!**
- Do not start any application; in case Remisol automatically starts, stop it,
- Back up the Remisol database. To do so, click on Windows **Start** button in the desktop lower left corner > **(All Programs if applicable) > 3 – Save Remisol Advance database**; the Remisol database will be backed up to its default location: **D:\MSSQL\Backup** under the name **Custom.dat**,
- Plug in the USB stick,
- Explore the USB stick folders and look for the one in which you have downloaded KB/Service Pack and the executable.
- Double click on the file SYSADM.EXE, the Windows session will close.
- Log on the machine again with the login with the highest administrative rights in your laboratory (sysadm)!**
- Double click on the Service Pack 3 installer located in the USB stick and follow onscreen instructions.
- Restart the Remisol client,
- Log on the machine again with the login with the highest administrative rights in your laboratory (sysadm)!**
- Double click on the patch / KB located in the USB stick to launch its installation – follow onscreen instructions,
- When completed, unplug the USB stick,
- Restart the client.

The client is now patched.



NORMAND-INFO

S.A.S.U. au capital de 17 500 001 €

7 rue Frédéric Degeorge – Bât B
2^{ème} étage - CS 10451
62028 ARRAS CEDEX
FRANCE

Tél : +33 (0)1.49.90.86.00
Fax : +33 (0)1.49.90.86.01

<http://www.normand-info.fr>

N° de SIRET 529 094 476 - Code APE 6201Z - RCS ARRAS 529 094 476 - N° Intra Communautaire FR 78529094476