


MicroScan Microbiology Systems



MAKING SPECIALTY A ROUTINE

Specialty identification panels for MicroScan microbiology systems

As yeast, fastidious organisms and anaerobes play a greater role in infectious disease, MicroScan specialty identification panels meet the challenge. Specifically, these panels:

- > Provide definitive identification in only four hours
- > Perform and interpret tests both manually and with a variety of MicroScan instruments
- > Deliver a cost-effective solution for supplemental testing
 - Simple setup
 - Easy-to-read Substrates
 - Automated interpretation through LabPro Information Manager or online by the Biotype Lookup tool

Rapid Yeast ID B1017-70

<i>Blastoschizomyces capitatus</i>	<i>Candida pintolopesii</i>	<i>Hansenula polymorpha</i>
<i>Candida albicans</i>	<i>Candida rugosa</i>	<i>Kluyveromyces lactis</i>
<i>Candida catenulata</i>	<i>Candida tropicalis</i>	<i>Pichia anomala</i>
<i>Candida famata</i>	<i>Candida viswanathii</i>	<i>Pichia farinosa</i>
<i>Candida glabrata</i>	<i>Candida zeylanoides</i>	<i>Prototheca</i> species
<i>Candida guilliermondii</i>	<i>Cryptococcus albidus</i>	<i>Prototheca wickerhamii</i>
<i>Candida inconspicua</i>	<i>Cryptococcus humicola</i>	<i>Rhodotorula glutinis</i>
<i>Candida kefyr</i>	<i>Cryptococcus laurentii</i>	<i>Rhodotorula minuta</i>
<i>Candida krusei</i>	<i>Cryptococcus melibiosum</i>	<i>Rhodotorula mucilaginosa</i>
<i>Candida lambica</i>	<i>Cryptococcus neoformans/gattii</i> complex	<i>Saccharomyces cerevisiae</i>
<i>Candida lipolytica</i>	<i>Cryptococcus terreus</i>	<i>Sporobolomyces salmonicolor</i>
<i>Candida lusitanae</i>	<i>Cryptococcus uniguttulatus</i>	<i>Trichosporon beigelii</i> complex
<i>Candida parapsilosis</i> group	<i>Geotrichum</i> species	

Rapid Haemophilus/Neisseria ID (HNID) Panel B1012-10B

<i>Aggregatibacter aphrophilus</i>
<i>Gardnerella vaginalis</i>
<i>Haemophilus haemolyticus</i>
<i>Haemophilus influenzae</i> Biotypes I, II, III, IV, V, VI, VII
<i>Haemophilus parainfluenzae</i> Biotypes I, II, III, IV
<i>Moraxella catarrhalis</i>
<i>Neisseria gonorrhoeae</i>
<i>Neisseria lactamica</i>
<i>Neisseria meningitidis</i>
<i>Neisseria mucosa</i>
<i>Neisseria</i> species



Rapid chromogenic panel databases

Rapid Anaerobe ID Panel B1017-2

Anaerobic cocci	Clostridia	Anaerobic	Anaerobic Non-Sporeforming Gram-Positive Bacilli (NSF-GPB)
<i>Acidaminococcus fermentans</i>	<i>Clostridium baratii</i>	<i>Bacteroides eggerthii</i>	<i>Actinomyces israelii</i>
<i>Anaerococcus prevotii</i>	<i>Clostridium bifermentans</i>	<i>Bacteroides fragilis</i>	<i>Actinomyces odontolyticus</i>
<i>Anaerococcus tetradius</i>	<i>Clostridium butyricum</i>	<i>Bacteroides ovatus</i>	<i>Actinomyces viscosus</i>
<i>Fingoldia magna</i>	<i>Clostridium cadaveris</i>	<i>Bacteroides thetaiotaomicron</i>	<i>Bifidobacterium dentium</i>
<i>Peptostreptococcus anaerobius/stomatus</i>	<i>Clostridium clostridioforme</i>	<i>Bacteroides uniformis</i>	<i>Eggerthella lentum</i>
<i>Peptoniphilus asaccharolyticus</i>	<i>Clostridium difficile</i>	<i>Bacteroides ureolyticus</i>	<i>Eubacterium limosum</i>
<i>Staphylococcus saccharolyticus</i>	<i>Clostridium histolyticum</i>	<i>Bacteroides vulgatus</i>	<i>Lactobacillus</i> species
<i>Veillonella parvula</i>	<i>Clostridium innocuum</i>	<i>Capnocytophaga</i> species	<i>Propionibacterium acnes</i>
	<i>Clostridium perfringens</i>	<i>Fusobacterium mortiferum</i>	<i>Propionibacterium granulosum</i>
	<i>Clostridium ramosum</i>	<i>Fusobacterium necrophorum</i>	
	<i>Clostridium septicum</i>	<i>Fusobacterium nucleatum</i>	
	<i>Clostridium sordellii</i>	<i>Fusobacterium varium</i>	
	<i>Clostridium sporogenes</i>	<i>Parabacteroides distasonis</i>	
	<i>Clostridium subterminale</i>	<i>Prevotella bivia</i>	
	<i>Clostridium tertium</i>	<i>Prevotella buccae</i>	
	<i>Clostridium tetani</i>	<i>Prevotella corporis</i>	
		<i>Prevotella disiens</i>	
		<i>Prevotella melaninogenica</i>	
		<i>Prevotella oralis</i> group	
		<i>Porphyromonas asaccharolytica</i>	
		<i>Porphyromonas gingivalis</i>	

Ancillary Products	Catalog Number	Yeast	Anaerobe	HNID
0.05N Sodium Hydroxide	B1015-3	•		
0.5% N,N-Dimethyl-Alpha-Naphthylamine	B1010-45A		•	•
0.8% Sulfanilic Acid	B1010-44A		•	•
Cover Tray	B1010-56B	•	•	•
Ehrlich's Reagent	B1015-6		•	
HNID® Indole Reagent	B1015-15			•
HNID® Inoculum Broth	B1015-16			•
Inoculum Water	B1015-2	•	•	
Mineral Oil, 60 mL	B1010-40		•	•
Peptidase Reagent	B1012-30B	•	•	•
Xylene	B1015-5		•	
Yeast Turbidity Standard	B1015-18	•		

For more information about specialty identification panels, please visit us at:
www.beckmancoulter.com/microbiology

©2017 Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo and the Beckman Coulter product and service names mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For Beckman Coulter's worldwide office locations and phone numbers, please visit www.beckmancoulter.com/contact

DS-52090

