



New Immunoassay Access Product Label Examples

Reagent Kit Label:

ACCESS HYBRITECH PSA			
<p>INFORMATION FOR USA ONLY Rx Only. 2 Reagent Packs R1, for the quantitative determination of PSA levels in human serum using the Access Immunoassay Systems. Reactive Ingredients: Paramagnetic particles coated with mouse monoclonal anti-PSA, BSA, 3.37 mL (R1a), Mouse monoclonal anti-PSA alkaline phosphatase (bovine) conjugate, BSA, protein (mouse), 3.1 mL (R1b).</p>			
CONTENTS 3.37 mL 3.1 mL	R1a R1b		

Calibrator Kit Label:

ACCESS HYBRITECH PSA CALIBRATORS			
<p>INFORMATION FOR USA ONLY Rx Only. For the quantitative determination of PSA levels in human serum using the Access Immunoassay Systems. Reactive Ingredients: BSA matrix 1 x 2.5 mL (S0), BSA, Human PSA 5 x 2.5 mL (S1-S5).</p>			
CONTENTS 1 x 2.5 mL 5 x 2.5 mL 2	S0 S1-S5 		

Control Kit Label:

ACCESS HYBRITECH PSA QC			
<p>INFORMATION FOR USA ONLY Rx Only. For monitoring system performance of the Access Hybritech PSA assay. Reactive Ingredients: BSA, Human PSA 3 x 5.0 mL (QC1-QC3).</p>			
CONTENTS 1 x 5.0 mL 1 x 5.0 mL 1 x 5.0 mL 2	QC1 QC2 QC3 		

Vial Label:

ACCESS Hybritech PSA CALIBRATOR		
	S0	
Hyb 0 ng/mL WHO 0 ng/mL		

Warning label (back of box)

